


Tagging Plasmid: pRPa^{TAG}



Expression cassettes

GFP_X / Myc_X

AAGCTT **ATG-GFP or Myc** TCTAGA **insert STOP** GGATCC
HindIII XbaI BamHI

X^{GFP}

TTAATTAAGCTT **ATG-insert** TCTAGA **GFP-STOP** GGATCC
PacI HindIII XbaI BamHI

X^{6Myc}

AAGCTTAATTA **ATG** TCTAGA **6Myc** CCTAGG **STOP** GGATCC
HindIII PacI XbaI AvrII BamHI

Promoter region

RRNA promoter +1 Tet operator KpnI Polypyrimidine tract
 ...AAGTAGCGCTTACGGCGTACCTATCAATGATAGAGTGGTACCCTGCACGCGCCTTCGAGTTTTTTTTCCTTTTCCCCATTTTTCAACTTGA
 AGA ACTTCAATTACACCAAAAAGTAAAATTCACAAGCTT -- Expression cassette
 SA HindIII

pRPa^{TAG}

For inducible expression of N/C-terminal human cMYC- (EQKLISEEDL) or N/C-terminal enhanced GFP-tagged proteins in *T. brucei* from a tetracycline-responsive RNA promoter.

- Derived from pLEW100 (Wirtz et al, 1999) and p2T7^{TAbue} (Alibu et al., 2005).
- High fidelity polymerase recommended.
- Integrates at the tagged *RRNA* spacer in 2T1/TAG^{PAC} *T. brucei* following digestion with *Ascl* (Alsford et al 2005), giving a transformation efficiency of $\sim 2.5 \times 10^{-6}$ (Alsford & Horn, 2007).

Cloning

GFP/MYC_X

To ensure that your gene is in frame with the tag, place the second codon downstream of the *Xba*I site, i.e. TCTAGA:[codon 2]:[codon]_n: [stop]:GGATCC.

X^{GFP}

To ensure that your gene is in frame with the tag, place the last but one codon upstream of the *Xba*I site, i.e. AAGCTT:[start]:[codon]_{n-1}:TCTAGA.

X^{6MYC}X

N-terminal tagging: clone your gene without a start codon via *Avr*II/*Bam*HI digestion (or without start/stop codons via *Avr*II)
i.e. CCTAGG[codon 2]:[codon]_n: [stop]:GGATCC.

C-terminal tagging: clone your gene without a stop codon via *Hind*III(or *Pac*I)/*Xba*I digestion (without start/stop codons via *Xba*I)
i.e. AAGCTT:[start]:[codon]_n:TCTAGA.

Use *Hind*III / *Bam*HI if you don't want the tag.

There are alternatives if the gene contains *Xba*I, *Avr*II and/or *Bam*HI:

Plasmid	<i>Xba</i> I, <i>Avr</i> II	Insert	<i>Avr</i> II, <i>Nhe</i> I, <i>Spe</i> I, <i>Xba</i> I
	<i>Bam</i> HI		<i>Bgl</i> II

Key features

- Complete sequences available.
- Hygromycin for stable selection.
- All vectors allow inducible expression using tetracycline (or analogues).
- Inducible cassette is independent of selectable marker.
- Modular nature allows tag or other components to be exchanged.
- Compatible with *T. brucei* cell lines expressing TetR and containing a tagged *RRNA* spacer, e.g. 2T1/TAG^{PAC} (Alsford et al, 2005).

Upon integration into *T. brucei*, the construct replaces TAG^{PAC} and generates a functional *HYG*^R at the previously tagged *RRNA* spacer. The operator binds Tet-repressor in the absence of tetracycline so the inducible RNA promoter is activated and tagged protein is expressed when tetracycline (1 μ g ml⁻¹) is added to the medium.

Detection:

cMYC	Mouse anti-cMYC, 9E-10 (Santa Cruz; IFA / western blotting) Mouse anti-cMYC, 4A6 (Upstate Biotechnology; western blotting only, as binds the mitotic spindle in <i>T. brucei</i>)
eGFP	Rabbit anti-GFP, IgG fraction (Molecular Probes; IFA and western blotting)

Other questions/comments, contact David Horn (david.horn@lishtm.ac.uk).