



Restriction Map and Multiple Cloning Site (MCS) of pDsRed-Monomer-C1. Unique restriction sites are in bold.

Description

pDsRed-Monomer-C1 is a mammalian expression vector that encodes DsRed-Monomer (DsRed.M1), a monomeric mutant derived from the tetrameric *Discosoma sp.* red fluorescent protein DsRed (1). DsRed-Monomer contains forty-five amino acid substitutions (listed on page 2). When DsRed-Monomer is expressed in mammalian cell cultures, red fluorescent cells can be detected by either fluorescence microscopy or flow cytometry 12–16 hr after transfection (DsRed-Monomer excitation and emission maxima = 557 nm and 585 nm, respectively). The DsRed-Monomer coding sequence is human codon-optimized for high expression in mammalian cells (2).

The multiple cloning site (MCS) in pDsRed-Monomer-C1 is positioned between the DsRed-Monomer coding sequence and the SV40 polyadenylation signal (SV40 poly A). Genes cloned into the MCS will be expressed as fusions to the C-terminus of DsRed-Monomer if they are in the same reading frame as DsRed-Monomer and there are no intervening stop codons. A Kozak consensus translation initiation site upstream of DsRed-Monomer increases the translation efficiency in eukaryotic cells (3). SV40 poly A signals downstream of the MCS direct proper processing of the 3' end of mRNA transcripts. The vector backbone also contains an SV40 origin for replication in mammalian cells expressing the SV40 T-antigen, a pUC origin of replication for propagation in *E. coli*, and an f1 origin for single-stranded DNA production. A neomycin resistance cassette (Neo^r)—consisting of the SV40 early promoter, the neomycin/kanamycin resistance gene of Tn5, and polyadenylation signals from the Herpes simplex virus thymidine kinase (HSV TK) gene—allows stably transfected eukaryotic cells to be selected using G418. A bacterial promoter upstream of this cassette expresses kanamycin resistance in *E. coli*.

Use

pDsRed-Monomer-C1 can be used to construct fusions to the C-terminus of DsRed-Monomer. If a fusion construct retains the fluorescent properties of the native DsRed-Monomer protein, its expression can be monitored by flow cytometry and its localization *in vivo* can be determined by fluorescence microscopy. The target gene must be cloned into pDsRed-Monomer-C1 so that it is in frame with the DsRed-Monomer coding sequences, with no intervening in-frame stop codons. The recombinant DsRed-Monomer vector can be transfected into mammalian cells using any standard transfection method. If required, stable transformants can be selected using G418 (4). pDsRed-Monomer-C1 can also be used as a cotransfection marker; the unmodified vector will express DsRed-Monomer.

We recommend using the DsRed-Monomer-C sequencing primer (see the Sequencing primer location section below) to sequence genes cloned adjacent to the 3' end of the DsRed-Monomer coding region.

For Western blotting, the BD Living Colors™ DsRed Polyclonal Antibody (Cat. No. 632397) can be used to recognize the DsRed-Monomer protein. However, to generate optimal results it may be necessary to use a higher concentration of antibody than recommended on the DsRed Polyclonal Antibody Product Analysis Certificate.

Location of features

- Human cytomegalovirus (CMV) immediate early promoter: 1–589

Enhancer region: 59–465; TATA box: 554–560

Transcription start point: 583

C→G mutation to remove *Sac* I site: 569

- Human codon-optimized DsRed-Monomer gene

Kozak consensus translation initiation site: 606–616

Start codon (ATG): 613–615; Stop codon: 1366–1368

Amino acid substitutions (DsRed→DsRed-Monomer)

GCC→GAC (Ala-2 to Asp) mutation: 616–618

TCC→AAC (Ser-3 to Asn) mutation: 619–621

TCC→ACC (Ser-4 to Thr) mutation: 622–624

AAG→GAG (Lys-5 to Glu) mutation: 625–627

AAC→GAC (Asn-6 to Asp) mutation: 628–630

CGC→CAG (Arg-13 to Gln) mutation: 649–651

ACC→TCC (Thr-21 to Ser) mutation: 673–675

GAG→TAC (Glu-26 to Tyr) mutation: 688–690

CGC→AAG (Arg-36 to Lys) mutation: 718–720

CAC→ACC (His-41 to Thr) mutation: 733–735

AAC→CAG (Asn-42 to Gln) mutation: 736–738

GTG→GCC (Val-44 to Ala) mutation: 742–744

AAG→CAG (Lys-47 to Gln) mutation: 751–753

GTG→GCC (Val-71 to Ala) mutation: 823–825

AAG→ATG (Lys-83 to Met) mutation: 859–861

AAG→ACC (Lys-92 to Thr) mutation: 886–888

GTG→TCC (Val-96 to Ser) mutation: 898–900

ACC→GAG (Thr-106 to Glu) mutation: 928–930

ACC→CAG (Thr-108 to Gln) mutation: 934–936

TCC→ACC (Ser-117 to Thr) mutation: 961–963

ATC→AAG (Ile-125 to Lys) mutation: 985–987

TCC→GCC (Ser-131 to Ala) mutation: 1003–1005

ATG→GCC (Met-141 to Ala) mutation: 1033–1035

GCC→CCC (Ala-145 to Pro) mutation: 1045–1047

CGC→AAG (Arg-149 to Lys) mutation: 1057–1059

CGC→CAG (Arg-153 to Gln) mutation: 1069–1071

CAC→TCC (His-162 to Ser) mutation: 1096–1098

AAG→CAC (Lys-163 to His) mutation: 1099–1101

CTG→ACC (Leu-174 to Thr) mutation: 1132–1134

GTG→TGC (Val-175 to Cys) mutation: 1135–1137

GAG→GAC (Glu-176 to Asp) mutation: 1138–1140

TCC→ACC (Ser-179 to Thr) mutation: 1147–1149

ATC→GTG (Ile-180 to Val) mutation: 1150–1152

ATG→AAG (Met-182 to Lys) mutation: 1156–1158

TAC→AAC (Tyr-192 to Asn) mutation: 1186–1188

TAC→CAC (Tyr-193 to His) mutation: 1189–1191

TCC→AAC (Ser-203 to Asn) mutation: 1219–1221
 ATC→GTG (Ile-210 to Val) mutation: 1240–1242
 CGC→CAC (Arg-216 to His) mutation: 1258–1260
 ACC→GCC (Thr-217 to Ala) mutation: 1261–1263
 GGC→GCC (Gly-219 to Ala) mutation: 1267–1269
 CAC→TCC (His-222 to Ser) mutation: 1276–1278
 CTG→GGC (Leu-223 to Gly) mutation: 1279–1281
 TTC→TCC (Phe-224 to Ser) mutation: 1282–1284
 CTG→CAG (Leu-225 to Gln) mutation: 1285–1287

- Multiple Cloning Site: 1288–1368
- SV40 early mRNA polyadenylation signal
Polyadenylation signals: 1507–1512 & 1536–1541; mRNA 3' ends: 1545 & 1657
- f1 single-strand DNA origin: 1604–2059 (Packages the noncoding strand of DsRed-Monomer)
- Bacterial promoter for expression of Kan^r gene
–35 region: 2121–2126; –10 region: 2144–2149
Transcription start point: 2156
- SV40 origin of replication: 2400–2535
- SV40 early promoter
Enhancer (72-bp tandem repeats): 2233–2304 & 2305–2376
21-bp repeats: 2380–3000, 2401–2421 & 2423–2440
Early promoter element: 2456–2462
Major transcription start points: 2452, 2490, 2496 & 2501
- Kanamycin/neomycin resistance gene
Neomycin phosphotransferase coding sequences:
Start codon (ATG): 2585–2587; stop codon: 3377–3379
G→A mutation to remove *Pst* I site: 2767
C→A (Arg to Ser) mutation to remove *Bss*H II site: 3113
- Herpes simplex virus (HSV) thymidine kinase (TK) polyadenylation signal
Polyadenylation signals: 3615–3620 & 3628–3633
- pUC plasmid replication origin: 3964–4607

Sequencing primer location

- DsRed-Monomer-C sequencing primer (5'-AGCTGGACATCACCAACCACAACG-3'): 881–904
Note: The DsRed1-C Sequencing Primer (Cat. No. 632388) **cannot** be used as a sequencing primer for pDsRed-Monomer-C1.

Propagation in *E. coli*

- Suitable host strains: DH5 α , HB101, and other general purpose strains. Single-stranded DNA production requires a host containing an F plasmid such as JM109 or XL1-Blue.
- Selectable marker: plasmid confers resistance to kanamycin (50 μ g/ml) to *E. coli* hosts.
- *E. coli* replication origin: pUC
- Copy number: high
- Plasmid incompatibility group: pMB1/Col E1

Excitation and emission maxima of DsRed-Monomer

- Excitation maximum = 557 nm
- Emission maximum = 585 nm

References

1. Matz, M. V., *et al.* (1999) *Nature Biotech.* **17**:969–973.
2. Haas, J., *et al.* (1996) *Curr. Biol.* **6**:315–324.
3. Kozak, M. (1987) *Nucleic Acids Res.* **15**:8125–8148.
4. Gorman, C. (1985) In *DNA cloning: A Practical Approach*, Vol. II. Ed. D. M. Glover. (IRL Press, Oxford, U.K.), pp. 143–190.

Note: The attached sequence file has been compiled from information in the sequence databases, published literature, and other sources, together with partial sequences obtained by BD Biosciences Clontech. This vector has not been completely sequenced.

Notice to Purchaser

This product is intended to be used for research purposes only. It is not to be used for drug or diagnostic purposes nor is it intended for human use. BD Biosciences Clontech products may not be resold, modified for resale, or used to manufacture commercial products without written approval of BD Biosciences Clontech.

The CMV promoter is covered under U.S. Patent Nos. 5,168,062, and 5,385,839 assigned to the University of Iowa Research Foundation.

This product is the subject of pending U.S. and foreign patents.

Not-For-Profit-Entities: Orders may be placed in the normal manner by contacting your local representative or BD Biosciences Clontech Customer Service at either 800-662-2566 or 650-424-8222, extension 1. BD Biosciences Clontech grants not-for-profit research entities a worldwide, non-exclusive, royalty-free, limited license to use this product for non-commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product or its components to third parties. Any other use of this product will require a license from BD Biosciences Clontech. Please contact the Product Manager for Cell Biology at either 800-662-2566 or 650-424-8222, extension 7816; or by e-mail at licensing@clontech.com.

For-Profit entities that wish to use this product in non-commercial or commercial applications are required to obtain a license from BD Biosciences Clontech. For license information please contact the Product Manager for Cell Biology at either 800-662-2566 or 650-424-8222, extension 7816; or by e-mail at licensing@clontech.com.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2004 BD