

Labeled Donkey Anti–Rabbit IgG Antibodies

Table 1. Contents and storage information.

Material	Amount	Concentration	Storage	Stability
Fluorophore-labeled donkey anti-rabbit IgG (H+L) antibodies	0.5 mL	2 mg/mL solutions in 0.1 M sodium phosphate, 0.1 M NaCl, pH 7.5, 5 mM sodium azide	<ul style="list-style-type: none"> • 2–6°C • Protect from light • Avoid freeze-thaw cycles 	When stored as directed, products are stable for at least 3 months.
Degree of labeling: Typically 2–8 fluorophore molecules per IgG molecule; the exact degree of labeling is indicated on the product label.				
Approximate fluorescence excitation/emission maxima: See Table 2.				

Introduction

Fluorescent donkey anti–rabbit IgG antibodies (Table 2) from Invitrogen are prepared from affinity-purified antibodies that react with IgG heavy chains and all classes of immunoglobulin light chains from rabbit. The Alexa Fluor® dyes to which these antibodies are conjugated provide for extraordinarily bright antibody conjugates. The donkey anti–rabbit IgG antibodies show minimum crossreactivity to bovine, chicken, goat, guinea pig, hamster, horse, human, mouse, rat, and sheep serum proteins. The approximate fluorescence excitation and emission maxima for each of the conjugates are shown in Table 2.

In addition to the secondary antibodies described in this product manual, Invitrogen prepares fluorescent conjugates of many other species-specific anti-IgG antibodies, as well as conjugates of avidin, streptavidin, NeutrAvidin® biotin-binding protein, protein A, and protein G. For more information about these products, visit probes.invitrogen.com or contact Technical Support.

At the time of preparation, the products are certified to be free of unconjugated dyes and are tested in a cytological experiment to ensure low nonspecific staining.

Table 2. Labeled donkey anti-rabbit IgG antibodies.*

Catalog #	Label	Ex †	Em †
A10039	Alexa Fluor® 350	346	442
A21206	Alexa Fluor® 488	495	519
A10040	Alexa Fluor® 546	556	573
A31572	Alexa Fluor® 555	555	565
A10042	Alexa Fluor® 568	578	603
A21207	Alexa Fluor® 594	590	617
A31573	Alexa Fluor® 647	650	668
A10043	Alexa Fluor® 680	663	690

* Minimum crossreactivity to bovine, chicken, goat, guinea pig, hamster, horse, human, rabbit, mouse, rat, and sheep serum proteins. † Approximate fluorescence excitation (Ex) and emission (Em) maxima, in nm, for conjugates.

Guidelines for Use

We recommend centrifuging the protein conjugate solution briefly in a microcentrifuge before use and add only the supernatant to the experiment. This step eliminates any protein aggregates that may have formed during storage, thereby reducing nonspecific background staining.

Because staining protocols vary with application, the appropriate dilution of antibody should be determined empirically. For fluorophore-labeled antibodies, a final concentration of 1–10 µg/mL should be satisfactory for most immunohistochemical applications.¹

Reference

1. *Short Protocols in Molecular Biology, 2nd Edition*, F.M. Ausubel et al., Eds., John Wiley and Sons (1992) pp. 14-24–14-30.

Product List

Current prices may be obtained from our website or from our Customer Service Department.

Cat. no	Product Name	Unit Size
A10039	Alexa Fluor® 350 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A21206	Alexa Fluor® 488 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A10040	Alexa Fluor® 546 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A31572	Alexa Fluor® 555 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A10042	Alexa Fluor® 568 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A21207	Alexa Fluor® 594 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A31573	Alexa Fluor® 647 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL
A10043	Alexa Fluor® 680 donkey anti-rabbit IgG (H+L) *2 mg/mL*	0.5 mL

Visit www.invitrogen.com/antibody for reagents and tools for antibody research..

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probetech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

**For country-specific contact information,
visit www.invitrogen.com.**

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2007, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.