

**MOUSE ANTI-HUMAN TOPOISOMERASE II α (Ki-S1)
MONOCLONAL ANTIBODY**

CATALOG NUMBER:	MAB4197 (formerly Roche Catalog Number 1742353)	QUANTITY:	100 μ g
LOT NUMBER:		CONCENTRATION:	1 mg/mL
SPECIFICITY:	In immunoprecipitation and Western blot experiments the Ki-S1 antibody recognizes a major protein of 170 kD which was identified as the α isoform of topoisomerase II (1, 2). In addition, it has been shown that the Ki-S1 antibody recognizes a carboxyterminal α -isoenzyme specific epitope missing in topoisomerase II β (1). In immunohistochemistry the Ki-S1 antibody shows strong nuclear staining only in proliferating cells. The epitope recognized by the antibody is also detectable in paraffin-embedded tissue sections (3). Accordingly, it has been shown that the expression of topoisomerase II α is strongly restricted to proliferating cells (4). The topoisomerase II α antigen is preferentially expressed during G ₁ , S, G ₂ and M phase of the cell cycle, while resting, non-cycling cells (G ₀ phase) lack topoisomerase II α . In addition, constantly proliferating cells (e.g. cell lines) react positively to topoisomerase II α during the entire cell-cycle. This specificity of Ki-S1 antibody for proliferating cells might make it a useful tool for determination of the proliferative fraction in solid tumors such as mammary carcinomas (2, 5–8) and gliomas (9).		
IMMUNOGEN:	Nuclear protein preparation from the human cell line U937.50		
ISOTYPE:	IgG _{2a}		
CLONE NAME:	Ki-S1		
APPLICATIONS:	Western blot: 1-10 μ g/mL Immunoprecipitation Immunocytochemistry: 5-10 μ g/mL Immunohistochemistry: 5-10 μ g/mL Flow cytometry Optimal working dilutions must be determined by end user.		
SPECIES REACTIVITIES:	Human		
FORMAT:	Purified immunoglobulin.		
PRESENTATION:	Liquid in 0.02M Phosphate buffer, 0.25M NaCl, pH 7.6 with 0.1% sodium azide.		
STORAGE/HANDLING:	Maintain at 2-8°C for up to six months from date of receipt.		

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

For research use only; not for use as a diagnostic.

© 2003-2005: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com